



Surgical sets

PANEPACK

... RIGHT DIRECTION WITH SMART SETS

 **Medana**
Medsurgical & Anaesthesia Devices

This top custom surgery set group includes not only drapes, gowns, dressings and accessory material, but also another tools and instruments used during the surgery procedures. These sets are produced solely as customized products. This set group makes a complex system and due to the optimal composition provide maximum utility value. The advantages of this solution are in outstanding time savings during the surgery intervention preparation, simplified supply, surgery process and waste disposal costs reduction.

UNIVERSAL SET

description	spec.	size	pcs/set
Adhesive drape	self adhesive	80x100 cm	2
Adhesive drape	self adhesive	160x240 cm	1
Adhesive drape	self adhesive	160x180 cm	1
Drape 2-ply		100x150 cm	1
OP-tape		10x50 cm	1
Table cover		150x260 cm	1
Mayo stand cover		80x140 cm	2
Surgical gown Standard	size L		1
Surgical gown Standard	size XL		3
Laparotomy sponges	X-RAY	45x45 cm	10
Gauze balls	X-RAY	30x30 cm	40
Gauze balls	X-RAY	9x9 cm	20
Gaze swabs	X-RAY	7x90 cm	10
Electrocautery pencil			1
Suction set Ch. 22	Ø 7,3 mm	length 210 cm	1
Scalpel blade #.21			1
Syringe 20 ml			1
Needle counter			1
Latex gloves, powder free	size L		2 pair



KNEE SET

description	spec.	size	pcs/set
Arthroscopy drape with elastic fenestration	Fen. 6x8 cm	200x320 cm	1
Drape 2-ply		100x120 cm	1
Leggings		30x60 cm	1
OP-tape		10x50 cm	2
Table cover		150x190 cm	1
Mayo stand cover	reinforced, yellow	80x140 cm	1
Surgical gown Comfort	size XL		2
Surgical gown Standard Plus	size L		1
Surgical gown Standard Plus	size XL		1
Gauze balls		40x40 cm	10
Laparotomy sponges	X-RAY	45x45 cm	10
Gauze swabs		10x10 cm	10
Gauze swab filled with cellulose wadding		15x160 cm	1
Knitted bandage hydrophilic		16cm x 10m	1
Wound dressing		10x20 cm	1
Suction set Ch. 30	∅ 8 mm	length 350 cm	1
Electrocautery pencil			1
Scalpel blade #20			3
Redon drain Ch. 12			1
Irrigation syringe	100 ml		1
Incise film		35x44 cm	1
Electrocautery tip cleaner			1
Wrapping sheet			1



HIP SET

description	spec.	size	pcs/set
U-Adhesive drape	Fen. 9x100 cm	290x290 cm	1
Adhesive drape	self adhesive	160x240 cm	1
Drape 2-ply		75x90 cm	1
Leggings		35x120 cm	1
OP-tape		10x50 cm	4
Table cover		150x190 cm	1
Fluid collection pouch		30x40 cm	1
Mayo stand cover	blue, reinforced	80x140 cm	1
Surgical gown Standard	size L		1
Surgical gown Standard	size XXL		1
Surgical gown Standard Plus	size XXL		2
Gauze balls		30x30 cm	5
Laparotomy sponges	X-RAY, pre-washed	45x45 cm	20
Gauze swabs		10x20 cm	20
Suction set Ch. 22	Ø 7,3 mm	length 210 cm	1
Scalpel blade #.21			1
Scalpel blade #.24			3
Needle counter			1
Redon drain Ch. 10			1
Redon drain Ch. 12			1
Vacuum bottle	400 ml		2



C-SECTION SET

description	spec.	size	pcs/set
C-section drape with fluid collection bag	bag 80x100 cm	200x300 cm	1
Baby towel		70x80 cm	1
OP-tape		10x25 cm	1
Table cover		150x190 cm	1
Mayo stand cover		80x140 cm	1
Surgical gown Standard	size L		2
Surgical gown Standard Plus	size L		1
Surgical gown Standard Plus	size XL		1
Gauze balls	X-RAY	40x40 cm	30
Laparotomy sponges	X-RAY, pre-washed	45x45 cm	15
Wound dressing		9x25 cm	1
Electrocautery pencil			1
Scalpel blade #22			1
Needle 18G	pink		1
Syringe 2,5ml	Luer-Lock		1
Baby cap			1
Isolation gown with knitted cuffs			1
Electrocautery tip cleaner			1



ANGIOGRAPHY SET

description	spec.	size	pcs/set
Femoral / radial angiography drape	4 Fen	220x350 cm	1
Drape 2-ply		50x70 cm	1
Drape 2-ply		75x100 cm	1
High absorbent under pad		60x90 cm	1
Table cover		150x150 cm	1
PE-cover		100x110 cm	2
PE-cover		140x140 cm	1
Surgical gown Standard	size XXL		1
Gauze swabs		10x10 cm	40
Scalpel blade #11			1
Needle 18G	pink		1
Needle 22G	black		1
Syringe 2ml			2
Syringe 5ml			1
Syringes 10ml	Luer-Lock		2
Syringe 20ml			1
Guide wire	J3	0,035x150 cm	1
Manifold 2-port. right, OFF 550PSI with rotating adaptor			1
Pressure monitoring line		100 cm	1
Fluid administration set		180 cm	1
Needle counter			1
Bowl transparent	150 ml		1
Bowl blue	250 ml		1
Bowl blue	500 ml		1



OPHTHALMIC SET

description	spec.	size	pcs/set
Ophthalmic drape with two F/c pouches		160x240 cm	1
Table cover		100x150 cm	1
Armboard cover		20x50 cm	2
Phaco tray cover		50x80 cm	1
Surgical gown Standard	size L		2
Gauze balls		15x15 cm	5
Gauze swabs		5x5 cm	5
Gauze swabs		7,5x7,5 cm	10
Eye pad			1
Kocher forceps, plastic			1
Stab knife	15°		1
Slit knife	angled	2,5 mm	1
V.E Cannula 23G			1
Cannula Hydrodissection 27G			1
Needle 27G	gray	0,4x20 mm	1
Syringes 2,5ml	Luer-Lock		2
Syringes 5 ml	Luer-Lock		1
Bowl transparent	150 ml		1
Wrapping sheet		50x50 cm	1



CRANIOTOMY SET

description	spec.	size	pcs/set
U-split drape reinforced		150x200 cm	1
Adhesive drape 2-ply	adhesive	75x90 cm	2
Adhesive drape 2-ply	adhesive	150x180 cm	1
Table cover		150x190 cm	1
Drape SMS		100x150 cm	1
Fluid collection pouch, adhesive	1-chamber	30x50 cm	2
Velcro holder 22x2,5 cm	adhesive		1
Mayo stand cover		80x140 cm	1
Surgical gown Standard	size L		2
Gauze balls	X-RAY	9x9 cm	10
Gauze balls	X-RAY	20x20 cm	8
Gauze balls	2-ply	30x30 cm	5
Gauze swabs		10x10 cm	40
Gauze swabs		10x20 cm	20
Elastic bandage		10 cm x 5 m	1
Electrocautery pencil			1
Suction tube 30CH		300 cm	1
Kocher forceps, plastic			2
Scalpel blade #11			1
Scalpel blade #20			1
Needle 18G	pink		1
Needle 21G	green		1
Syringes 20ml			3
Fluid administration set		140 cm	1
Bowl blue	250 ml		2
Bowl blue	500 ml		1
Plastic tray blue		27x25 cm	1



JAW SET

description	spec.	size	pcs/set
Adhesive drape with Fen. 9x12 cm	adhesive	120x280 cm	1
Adhesive drape with Fen. 10x15 cm, 2-ply	adhesive	75x100 cm	1
Fluid collection pouch, adhesive	1-chamber	30x40 cm	1
Velcro holder 22x2,5 cm	adhesive		1
Table cover		100x150 cm	1
Camera drape		14x250 cm	1
Light handle cover universal			2
Mayo stand cover		80x140 cm	1
Drill cover		7x170 cm	1
Surgical gown Standard	size L		1
Surgical gown Standard	size XXL		2
Gauze balls	X-RAY	20x20 cm	4
Gauze balls		40x40 cm	3
Gauze swabs	X-RAY	7,5x7,5 cm	30
Gauze swabs	X-RAY	10x12,5 cm	40
Gauze swabs rolled	X-RAY	7x90cm	1
Suction set CH25, tapered tip	∅ 6 mm	length 350 cm	1
Syringes 2ml			4
Syringes 20ml			2
Extension tubing		1,8mm x 180 cm	2
Galipot transparent	125 ml		1



The logo for PANEP features the word "PANEP" in a bold, blue, sans-serif font. The text is centered within a red oval shape that is slightly open at the top and bottom.