



PAJUNK®

Rectus Sheath Kits


*Catheter System for a Continuous
Technique and Single Shot*

Rectus Sheath Kits

For Continuous Technique and Single Shot

Rectus sheath blocks have been described for many years. Abdominal field blocks since 1899 by Schleich¹ with the development of ultrasound techniques, rectus sheath catheters can now be placed with ultrasound guidance^{2,5}. Rectus sheath catheters have shown as part of an enhanced recovery/multimodal analgesia technique giving low pain scores and low opiate requirements.^{3,4}

PAJUNK® offers continuous kits and single shot needles including TuohySono, the ultrasound visible needle.



TuohySono Needle

The TuohySono needle from Pajunk combines proven Tuohy tip with the unique reflection of the patented "Cornerstone Reflectors". The needles are now available in individual packaging.



1 Schleich CL. Schmerzlose Operation. Berlin: Springer, 1899:240-1.

2 Sandeman DJ, Dille AV. Ultrasound-guided rectus sheath block and catheter placement. ANZ Journal of Surgery 2008; 78: 621-3.

3 Cornish P, Deacon A. Rectus sheath catheters for continuous analgesia after upper abdominal surgery. ANZ Journal of Surgery 2007; 77: 84.

4 Webster, K, Hubble, S. Rectus sheath analgesia in intensive care patients: technique description and case series. Anaesthesia and Intensive Care 2009; 37: 855.

5 Webster, K. Ultrasound guided rectus sheath block – analgesia for abdominal surgery. Update in Anaesthesia 2010; 26; 126 Peng Ultrasound-Guided Interventional... 2013, 38 (4): 27.

Rectus Sheath

Product	Item No. LUER	Item No. NRFit	PU
Rectus Sheath Kit	0731153-49	0731163-49	10
1 x Touhy needle 16G x 80mm plastic stylet 2 x catheters 19G x 50cm with closed tip 2 x Clamping Adaptors (including 2 x White Luer Caps) 1 x Y connector (including 3 x White Luer Hub Caps connected to each end of the tubing) 1 x filter 2 x (Rectus Sheath) labels –Yellow- 2 x FixoCath 1 x FixoLong 1 x Luer/NRFit Cap White			



LUER





NRFit

TuohySono

Product	Size	Item No.	PU
TuohySono			
with Cornerstone Reflectors, metal stylet and detachable plastic wings	22G x 120mm	1085-4E120	10
	22G x 100mm	1085-4E103	25
	22G x 90mm	1085-4E090	25
	22G x 70mm	1085-4E070	25
	22G x 50mm	1085-4E050	25
	21G x 120mm	1085-4F120	10
	21G x 90mm	1085-4F090	25
	21G x 50mm	1085-4F050	25
	20G x 100mm	1085-4G103	25
	18G x 150mm	1085-4K150	10
	18G x 120mm	1085-4K120	10
	18G x 90mm	1085-4K090	25
	18G x 50mm	1085-4K050	25
	17G x 150mm	1085-4M150	10
17G x 90mm	1085-4M090	25	

Get the app!



PAJUNK® GmbH
Medizintechnologie
Karl-Hall-Strasse 1
D-78187 Geisingen / Germany
Phone +49 (0) 77 04/92 91-0
Telefax +49 (0) 77 04/92 91-6 00
www.pajunk.com

PAJUNK® UK Medical Products Ltd
Unit C1 The Waterfront,
Goldcrest Way, Newburn Riverside,
Newcastle upon Tyne, NE 15 8NY,
England
Phone +44 191 264 7333
pajunk.co.uk