

Rythmic™ Evolution

PRODUCT
INFORMATION











Reference	<p>KP5.04.125.4 Rythmic™ Evolution (blue)</p> <p>KP5.04.136.4 Rythmic™ Evolution (blue rechargeable)</p> <p>KP5.04.261.3 Mini Rythmic™ Evolution Lockable (blue)</p> <p>KP5.04.265.3 Mini Rythmic™ Evolution Lockable (blue rechargeable)</p> <p>KP5.04.141.4 Rythmic™ Evolution Organiser 500 (blue)</p>	<p>KP5.04.144.5 Rythmic™ Evolution Organiser 501 (blue rechargeable)</p> <p>KP5.04.269.1 Mini Rythmic™ Evolution Lockable (blue rechargeable) UK</p> <p>KP5.04.190.1 Rythmic™ Evolution (blue rechargeable) UK</p> <p>KP5.04.151.5 Rythmic™ Evolution Organiser 501 (blue rechargeable) UK</p> <p>KP5.04.186.1 Rythmic™ Evolution (blue rechargeable, No key) UK</p>
Delivery modes	<ul style="list-style-type: none"> - Continuous - Bolus only - Continuous + bolus - Programmed Automatic Bolus (PAB) with combined PCA 	
Flow rate and Accuracy	0.1 to 100 ml/hr (configurable maximum rate) at +/- 5%	
Programming units Concentration	ml/mg/µg 0.1 to 99.9 mg/ml and 0.1 to 999.9 µg/ml	
Volume to infuse	1 to 9999 ml	
Bolus dose	0 or from 0.1 to 60 ml with 0.1 ml increment	
Bolus flow rate	100 ml/hr, adjustable from 50 to 100 ml/hr with 10 ml/hr increment	
Lockout interval	0 to 999 minutes	
Special functions	<p>Protocol libraries (60 protocols per library) - Dose limit / x hours and Bolus limit / x hours - Loading dose - Clinician bolus override - Mandatory priming (enabled, warning only, disabled) - Automatic priming through the pump - Lockable pump keyboard</p>	
Therapy monitoring	Recording of patient pain scores with visual analogue scale (VAS) - VAS available in the pump history - Printable and computer downloadable - Ready to Connect to www.micrelcare.net for therapy monitoring through the web	
KVO	0.1 to 2.0 ml/hr	
Automatic restart after occlusion release	Yes	
Safety features	Air in line elimination filter - Free flow protection - Programming within protocol limits - Lockable protective cover for drug container - Code protected stop of infusion	
Alarms & warnings	Occlusion (low 0.8 bar, middle 1.2 bar, high 1.6 bar) - Near end and end of infusion - Low and depleted battery - Dose limit, bolus limit - Open mechanism - Preventive maintenance warning	
Access code levels	3	
Programming	Protocol libraries programmed with 'Rythmic™ Evolution Therapy Manager' computer software Direct programming on the pump possible	
Languages	Bulgarian, Croatian, Czech, Danish, Dutch, English, Finnish, French, German, Greek, Hungarian, Italian, Lithuanian, Norwegian, Polish, Portuguese, Russian, Serbian, Slovenian, Spanish, Swedish	
History, infusion graph / printout and PC download	Yes, including 6500 time stamped events event log - 'Plug and play' print facility	
Power supply	100 to 240 V, 50-60 Hz (with optional adaptor KS5.01.615.4)	
Preventive maintenance	Every 3 years	



Includes bolus, battery, carrying bag & keys
(carrying bag applies for Rythmic™ Evolution and Mini Rythmic™ Evolution models)

 micrel

DEHP free set range

Rythmic™ Administration Set range for Rythmic Evolution Systems			
Mini Full Set 100 mL KE1.EE.177.0		Spike Set S KM1.EE.155.1	
Full Set 200 mL KE1.EE.185.1		Spike Set S Y KM1.EE.187.2	
Y Mini Full Set 100 mL KE1.EE.195.1		Short Set 100 mL KE1.EE.158.0	
Y Full Set 200 mL KE1.EE.196.0		Luer Set KE1.EE.197.0	
Spike Set KM1.EE.099.2		Valve Extension Set with filter 0,2 um KE1.EE.220.1	

Accessories for Rythmic™ Evolution Systems					
	Rythmic™ carrying bag Included in all Rythmic™ Evolution systems KM1.YY.415.2		Carrying bag 500 ml for Rythmic™ and IP Connect KM1.YY.318.2		Mini Rythmic™ Bag 500 ml KM1.YY.365.1
	Mini Rythmic™ Bag Included in all Rythmic™ Evolution systems KM1.YY.400.1		Rythmic™ XS Carrying Shoulder Bag single use KM1.YY.566.1		Rythmic™ XS Carrying Waist Bag 250 mL KM1.YY.506.2
	IP Connect Pack KP5.04.259.2		Rythmic™ Evolution Bolus Handset KS5.04.505.2		Rythmic™ Power Adaptor plug type C (EU) (Available also UK and USA Nema) KS5.01.615.4
	Mini Rythmic™ Pole Clamp KS5.04.017.5		Rythmic™ holder with Pole Clamp KS5.01.430.3		Rythmic™ Organizer Pole Clamp KS5.01.524.3

96081 en EN Rev.6 Rel Date 2023.11.21