

Technical Specifications

Flow	High speed: 1200ml/min
	Low speed: 500ml/min
Power	8x1.5V cell
Size & Weight	Extremity bag: 700x430 mm
	Standard trunk bag: 300x350 mm
	XL trunk bag: 400x430 mm
	Weight- irrigator: 165g
Accessories	Standard tip
	Flexible tunnel tip
	3 into 1 fluid adaptor
	100g solidifier to solidify approx.
	3L normal saline
Operating conditions	Temperature: 10°C - 40°C
	Relative Humidity: 0% - 90%
Shipping and storage conditions	Temperature: -20°C - 49°C
	Relative Humidity: 0% - 90%

Part Number	Description	UOM
Kit		
EK-100K	Irrigator, with 3 into 1 fluid adaptor, standard Tip and Extremity Bag (No Granules) ALL items Sterile in single blister pack	1/box
AB-100K	Irrigator, with 3 into 1 fluid adaptor, standard Tip and Trunk Bag(No Granules) standard Size, ALL items Sterile in single blister pack	1/box
AB-100K-XL	Irrigator, with 3 into 1 fluid adaptor, standard Tip and Trunk Bag (No Granules) Extra Large Bag Size, ALL items Sterile in single blister pack	1/box
PI-200H	Irrigator, with 3 into 1 fluid adaptor, ALL items Sterile in single blister pack	1/box
Bag and Tip		
EB-100	Un-Sterile Extremity Bag WITH Granules, and Sterile Tip	1/box
AB-100	Un-Sterile Trunk Bag size 300mm x 350mm WITH Granules, and Sterile Tip	1/box
AB-100-XL	Un-Sterile Trunk Bag Extra Large, bag size 400mm x 430mm WITH Granules, and Sterile Tip, 1 piece per box	1/box
EB-100-3	Un-Sterile Extremity Bag WITH Granules, and Sterile Tip	3/box
AB-100-3	Un-Sterile Trunk Bag size 300mm x 350mm WITH Granules, and Sterile Tip	3/box
AB-100-3-XL	Un-Sterile Trunk Bag Extra Large, Bag size 400mm x 430mm WITH Granules, and Sterile Tip	3/box
EB-100-10	Un-Sterile Extremity Bag WITH Granules and Sterile Tip	10/box
AB-100-10	Unsterile Trunk Bag, size 300mm x 350mm WITH Granules and Sterile Tip	10/box
AB-100-10-XL	Unsterile Trunk Bag Extra Large, Bag size 400mm x 430mm WITH Granules and Sterile Tip	10/box
Tip and Adaptor		
FT100-10	Sterile Flexible Tip	10/box
ST100-10	Sterile Standard Tips	10/box
HY100-10	Sterile Hydra Adaptor	10/box

For more information, please refer to the instruction manual.



Wound Care Solutions Limited

Matrix Business Centre, Nobel Way, Dinnington, Sheffield, South Yorkshire, S25 3QB, United Kingdom Tel: +44(0) 1909 547232 Fax:+44(0) 1909 547233 www.woundcs.com